



# Illumina Sample Prep – Best Practices and Troubleshooting

Rey Sequerra, Field Applications  
Scientist

© 2010 Illumina, Inc. All rights reserved.

Illumina, illuminaDx, Solexa, Making Sense Out of Life, Oligator, Sentrix, GoldenGate, GoldenGate Indexing, DASL, BeadArray, Array of Arrays, Infinium, BeadXpress, VeraCode, IntelliHyb, iSelect, CSpPro, GenomeStudio, Genetic Energy, HiSeq, and HiScan are registered trademarks or trademarks of Illumina, Inc. All other brands and names contained herein are the property of their respective owners.


**illumina**<sup>®</sup>

# Outline

- ▶ TruSEQ Best practices
  - Starting material
  - DNA fragmentation
  - Reagent Handling
  - Agarose gel size selection
- ▶ Bioanalyzer Trace as a Troubleshooting tool
  - Commonly seen problems
- ▶ TruSEQ Small RNA Troubleshooting
- ▶ Nextera

# Library Preparation

- ▶ Prepares sample nucleic acid for sequencing
  - Fragmenting
  - Generates double-stranded DNA (if necessary)
  - Flanks with Illumina adapters
- ▶ All preparation ends with the same product
  - Double-stranded DNA with insert to be sequenced flanked by adapters
- ▶ Same protocols all both platforms




Illumina's Library Preparation Workflow

# Overview of TruSeq RNA Workflow

*Master-mixed reagents, optimized adapters, and a flexible workflow for simple method to prepare RNA for NGS*


- ▶ Isolate mRNA from Total RNA with poly-A selection
- ▶ Randomly Fragment RNA
- ▶ Random Prime mRNA → cDNA
- ▶ Make 2nd Strand cDNA
- ▶ End repair, phosphorylate, A-tailing
- ▶ Ligate sequencing adapters (w/ indexes)
- ▶ Gel free size selection
- ▶ Enrich with 15 cycles of PCR
- ▶ Grow Clusters and Sequence



# TruSeq DNA & RNA v2 Sample Prep

*Additional 12 indexes available for a total of 24!*

- ▶ Each TruSeq DNA and RNA v2 kits:
  - No changes to workflow or price, increase in index capabilities
  - Fill volumes and new consumables to support automation
  - Each kit contains 12 of 24 unique indexes
  - Each index reaction sufficient for 8 individual samples



Adapter

# TruSeq Exome and Custom Enrichment

*Simple, integrated and cost-effective*



- ▶ Simple and scalable workflow
  - Plate-based processing and master-mixed reagents
  - Gel-free; integrated with TruSeq DNA Sample Prep
  - 1ug DNA input
  - 3 days total time with less than 4hr hands-on
- ▶ Pre-enrichment sample pooling
  - Dramatically reduces FTE time and increases throughput

	Exome	Custom
Region size	62Mb	500kb-25Mb
Pre-enrichment pooling	6 spls	12 spls
Read length	2 x 100	2 x 50-75
Min order	48 spls	48 spls
Demo/Test kit available	Exome	Cancer panel



# Best Practices: Starting Material

- ▶ Use recommended input amounts
- ▶ Quantify with Qubit, Picogreen, or Ribogreen
- ▶ Assess sample quality
  - 260/280 ratio and gel electrophoresis for DNA
  - RIN score for RNA (Bioanalyzer)
- ▶ Determine sample source for RNA protocols
  - Total RNA, purified mRNA
  - Origin species: no poly-A tail, dicer like processing





Figure 8 Analysis of Genomic DNA Sample Integrity



# Best Practices:

## Input Recommendations

	Recommended		Recommended
gDNA (SR;PE; Indexed)	500ng - 5 µg*	TruSeq DNA	1ug genomic
Nextera	50ng		
Targeted Resequencing	1 - 5µg	Enrichment	500ng DNA library
Mate Pair gDNA	10 - 20µg		
ChIP-Seq	10 - 50 ng		
mRNA	0.1 - 4 µg total RNA	TruSeq RNA	1ug total RNA
Small RNA	1 µg total RNA	TruSeq small RNA	1ug total RNA

\*100 ng has been shown to work



# Best Practices:

## Starting Material, Enrichment

- ▶ Assess library quality and size
  - Run agarose gel
  - Run Bioanalyzer
- ▶ Determine library quantity
  - Fluorometric assay, like Qubit dsDNA Assay
  - Use 2ul DNA library in 200ul Qubit buffer
- ▶ Use 500ng of purified DNA library in Enrichment prep


# Best Practices:

## Fragmentation

- ▶ Fragmentation method: Covaris, nebulizer, sonicator
  - Were adjustments made for sample loss
  - Recommended to run a bioanalyzer trace, especially the first couple of times using a new fragmentation protocol
- ▶ Covaris fragmentation can be impaired by excessive air bubbles or airgap
  - Spin tube down before fragmentation to minimize risk
- ▶ Can sonicate up to 10ug depending on protocol
  - Ensure that only 1ug goes into TruSeq protocol (Follow recommendations on protocol input material)

# Best Practices: Fragmentation


- ▶ Example of good Covaris fragmentation



# Fragmentation Problems

## DNA Tailing


- ▶ Incomplete shearing caused by too much DNA input



- ▶ Adhere to recommended starting amounts

## Increased size range



- ▶ Covaris problem



- ▶ Check Covaris with DNA of known quality (like lambda)

# Fragmentation Problems

- ▶ Uneven shearing
- ▶ Poor DNA quality
- ▶ DNA eluted in wrong buffer (not TE)
- ▶ Covaris issues
  - Volume too low. Allows air gap and causes inconsistent fragmentation
  - Water level is too high or low
  - Temperature of water bath not between 6-8°C
  - Insufficient degassing



# Best Practices:

## Master Mix Reagent Handling

- ▶ Minimize freeze-thaw cycles
  - Dispense the reagent into aliquot
  - Beware: more aliquot -> reduced excessive reaction
- ▶ Reagents must be added in ORDER indicated in the TruSeq Sample Prep Guide
- ▶ A-Tailing Mix (ATL) and Ligation Mix (LIG) are viscous



# Best Practices:

## AMPure XP Beads


- ▶ Clean up steps have only been verified using a 96-well PCR or MIDI plate. Comparable performance is not guaranteed when using a microcentrifuge tube or other formats
- ▶ Ensure beads are well suspended, not settled
- ▶ Carefully pipet up and down to mix to minimize potential bead loss
- ▶ Utilize freshly prepared 80% EtOH for the washes
- ▶ While drying, allow the plate to remain on the magnetic stand to prevent potential bead loss due to electrostatic forces




# Common Questions

## AMPure XP Beads

- ▶ Can I use columns in place of AMPure beads?
  - Not supported in the TruSeq protocol, but can work
  - Columns will not remove smaller DNA species, while the SPRI beads will remove everything below 150bp
  - Gel will not run as clean due to retention of the adapters
- ▶ Do the SPRI steps have to be done before the gel?
  - Gel size selection step will remove the adapters, but it is still recommended to do a SPRI cleanup before running the gel
  - This is because excess adapter dimers may cause the gel to run less reproducibly



# Agarose Gel for Size Selection



- ▶ Familiar
- ▶ Cheap
- ▶ Very tight selection
- ▶ Can take several slices
- ▶ Not readily automatable

## Tips:


- ▶ Take care loading the sample – to avoid losses and cross contamination
- ▶ Don't overload (keep to less than 3ug/lane)

# DNA Migration Through Gels


- ▶ To ensure accurate size selection, the protocol should be run exactly as listed in the user guide
- ▶ To minimize contamination, separate samples on the gel or run one sample per gel
- ▶ Ladders
  - We recommend a Promega ladder because it produces nice, sharp bands
    - BenchTop 100 bp DNA ladder, Promega, part # G829B
  - Do not overload your DNA ladder as this can lead to skewed migration patterns
- ▶ Pre- Vs. Post Gel Staining
  - We recommend a SyBr Gold Nucleic acid gel stain from Invitrogen
    - part # S11494

# Pre and Post Staining (SyBr Gold at 1x)

**Pre-electrophoresis Staining**




Selected 300 – 400 bp




Selected 400 – 500 bp

**Post-electrophoresis Staining**





Selected 300 – 400 bp



Selected 400 – 500 bp

# TruSeq DNA gel size selection

- ▶ Invitrogen E-gels have been shown to give inaccurate size selection
  - ▶ The selected product appears to be the right size after the gel but before amplification
  - ▶ After amplification, the product has a smaller apparent size and often shows a second larger band





# Troubleshooting:




## Size Selection Step

- ▶ Gather gel images before and after gel cutting from customer
- ▶ Ensure size range cut is what was intended for project
  - Make sure customer is running the correct ladder
- ▶ Determine size, type and staining method of gel
  - Gel should be 150ml, low-range ultra agarose prestained with SYBR gold
- ▶ Ensure that customer did not overload gel if starting with more input material than recommended




# Bioanalyzer trace as a troubleshooting tool

# Expected Bioanalyzer Traces



# Presence of peak <125bp: Adapter dimers



**Primer  
dimers**

**PCR artifact**

**DNA library (50% of total DNA)**

From: Agilent's Bioanalyzer Applications for Next-Gen Sequencing: Updates and Tips March 1, 2011

# Presence of peak <125bp: Adapter dimers

## Cause

- ▶ Error or omission of DNA gel cut
  - Proper cut should remove <150 bp
- ▶ Improper AMPure bead ratios
  - Proper ratios should remove <150bp
- ▶ Inefficient Ligation
  - Too much input DNA
  - Too little input DNA

## Effect

- ▶ Dimers will form clusters and be sequenced. If ratio is low compared to library, may not be a problem but some reads will be dimers.
- ▶ Cause low diversity sample with improper matrix/phasing info and low clusters passing filter if clustered at higher densities.

- Primers <100bp cannot cluster or be sequenced, but can bind to flowcell

# Libraries of Unexpected Size

## Cause

- ▶ Recommended gel conditions not used
- ▶ Check before and after gel purification pictures, if available



## Effect

- ▶ Other types of gels and staining methods can shift size ranges leading to incorrect band excision



# Libraries 2x larger than expected: PCR artifact

- ▶ Possible single stranded species
- ▶ Overamplification




From: Agilent's Bioanalyzer Applications for Next-Gen Sequencing: Updates and Tips March 1, 2011

# Additional Peaks:

## SPRI bead carryover

- ▶ SPRI bead carry over following PCR Clean Up step



- ▶ Use a strong magnet for bead separation
- ▶ Pipet carefully during elution to avoid disturbing bead pellet




# TruSeq Small RNA Troubleshooting

# TruSeq Small RNA Prep:

## Appropriate bands

- ▶ miRNA libraries 145-160bp
- ▶ 18bp insert = 136bp library
- ▶ 25bp insert = 143bp library
- ▶ All consistently run about 7bp higher on gel and Bioanalyzer


**Figure 12** Small RNA Library from Total RNA Samples



# TruSeq Small RNA Prep:

## Inappropriate bands

- ▶ <100bp: Normal, Primers or primer dimers
- ▶ 120-125bp: Amplified adapter dimer
  - Low input increases dimer
- ▶ 130-138bp: Artifact band, Adapter concatamer
  - Reduced by increasing temp of RT
  - Seen at 44-47', gone at 50'
  - Make sure ice->cycler->ice
- ▶ 144-150bp: PRODUCT. Cut band
- ▶ 150-155bp: pre-miRNA, piRNA
- ▶ 200bp smear: tRNA + adapters





# TruSeq Small RNA Prep

## Gel best practices

- ▶ Gels and gel running buffer should be fresh
  - Commercially available gels are recommended
- ▶ Don't run gel so fast that it warms up
  - Warmed up gels will warp
  - May need to optimize running conditions on individual apparatus
- ▶ Premix stain to proper concentration before adding to gel
  - Adding concentrated stain directly can cause blobs of stain that cant be washed away
- ▶ Don't overload gel
  - Load 25ul of post-PCR volume
  - More or less causes the band to shift up or down
- ▶ Don't mix by vortexing
  - Vortexing causes fuzzier bands and lower yields



# TruSeq Small RNA Prep

## Loss of yield

- ▶ No band between 145-160bp
  1. Mostly caused by formation of secondary structure in adapter oligos that prevents ligation
- ▶ Solution: Verify 70°C incubation of adapters was performed
- ▶ Band present, and excised, at 145-160bp but lost in subsequent steps
  1. Gel disruption or elution
  2. Ethanol precipitation (if done)
- ▶ Solutions:
  1. Practice precipitation with 25bp DNA ladder
  2. Use freshly diluted ethanol and fresh sodium acetate
  3. Final pellet doesn't adhere well to tube wall – do NOT lose during washes!



Nextera

# Nextera DNA Sample Prep

*Sequencing's fastest and easiest sample prep*


- ▶ Rapid and streamlined workflow
  - Complete protocol in < 90 min
  - Single well enzymatic reaction both fragments and adds adapter
  - no mechanical fragmentation/shearing required
  - Easy four step protocol
- ▶ Ultra-Low DNA input of only 50 ng
- ▶ High Throughput, High Indexing
  - Prepare up to 96 samples per batch
  - Master-mixed reagents and automation-friendly configurations
  - 96 indices supported on all Illumina sequencers



# Step 1: Tagmentation of template DNA

Transposomes


Genomic DNA



Tagmentation



## Step 2: PCR to add adapters and indices



# Step 3: Cleanup and Sequence

Tagmentation



Reduced-Cycle  
PCR Amplification



Sequencing-Ready Fragment

# Improvements to Illumina's Nextera DNA Preparation Kit compared to legacy kit from Epicentre


- ▶ Optimized PCR
  - Reagents now included in the kit
  - Less cycles
  - New polymerase
  - Reduced GC bias and error rates
- ▶ 96 indices supported on all ILMN sequencers
- ▶ Improved master mixes – easier protocol with less tubes
- ▶ Optimized tagmentation reaction
  - shift insert size distribution to support longer 2x150 reads on MiSeq
- ▶ Plate-based protocol to support larger batch sizes





# General Troubleshooting Strategy

- ❖ Typical complaint is about sample 'loss'
- ❖ Is the sample completely gone or are there unexpected peaks in final trace?
- ❖ If the sample is gone, likely cleanup problem
- ❖ If there are unexpected peaks, likely tagmentation or PCR problem
- ❖ Confirm exactly 50ng of high quality genomic DNA was used
- ❖ Confirm which step in the protocol is the problem, run diagnostic traces after each step




# Sample Input

## ❖ General Comments:

- ❖ We require 50ng of high quality genomic DNA of at least 300bp
- ❖ Sample concentration can be measured with Qubit and quality can be assessed on a gel or trace




## ❖ Failure Modes:

- ❖ over fragmentation --- too little sample
- ❖ under fragmentation --- too much sample
- ❖ poor performance in PCR --- degraded sample
- ❖ loss in coverage within 50bp of the ends of the template --- template too short






# Sample Input

- ❖ Failure Mode: over fragmentation --- too little sample




# Sample Input

- ❖ Failure Mode: under fragmentation --- too much sample



# Tagmentation

- ❖ General Comments:
  - ❖ This is a 5 min. incubation at 55C
  - ❖ The protocol lists a hold at 10C . . . This means bring the sample to 10C and proceed **immediately** to Zymo cleanup. The transposomes are active at 10C and will continue fragmenting the sample.
- ❖ Failure modes:
  - ❖ inhibition of transposase (under fragmentation) --- presence of ethanol in sample
  - ❖ overfragmentation --- incorrectly calibrated thermal cycler: too hot
  - ❖ overfragmentation --- incorrect incubation time: too long




# Zymo

## ❖ General Comments:

- ❖ This is a simple column purification.
- ❖ A bioanalyzer trace of 1uL undiluted eluate should show a size range between 150bp and 1kb



## ❖ Failure modes:

- ❖ low yield from PCR --- ethanol not added to wash buffer: transposomes may remain on DNA and inhibit PCR reaction
- ❖ low yield from clean- up --- spinning at wrong speed and/or using incorrect volumes




# Zymo

## ❖ Failure Mode: transposome inhibition of PCR



- ❖ The transposomes and both sets of PCR primers will add ~135bp to the template
- ❖ It is unlikely you will ever see full adapter+primer dimers
- ❖ PCR primer dimers will likely run at ~55-60bp






# PCR

## ❖ General Comments:

- ❖ This is the most complicated stage of the protocol: a 4 primer, 5 cycle PCR reaction
- ❖ A very common failure mode is forgetting to add the two index primers
  - ❖ Two primers are in the PPC (PCR Primer Cocktail), the other two are in the index kit (purchased separately)

## ❖ Failure modes:


- ❖ low yield, contamination --- pipetting inaccuracies
- ❖ low yield --- forgetting to spin before PCR
- ❖ low yield, dimer formation (55-60bp) --- wrong PCR setting on thermal cycler
- ❖ low yield, dimer formation (55-60bp) --- miscalibration of thermal cycler
- ❖ low yield, evaporation of sample --- not heating lid on thermal cycler



# Ampure

- ❖ Important to use freshly prepared 80% ethanol, old ethanol will decrease wash efficiency
- ❖ As always, it is very important to let the supernatant clear before moving on
- ❖ Incorrect usage of magnetic stand can lead to sample loss and bead carryover
- ❖ Incorrect bead handling can cause contamination and/ or sample loss

Example of bead carryover



Sample Input

Tagmentation

Zymo

PCR

Ampure



Questions?